



BioGears Models: Multi-Trauma Physiology

A. Baird¹, J. Carter¹, L. Marin¹, M. McDaniel¹, N.
Tatum¹, S. White¹

1. Applied Research Associates Inc.

Outline

- Multi-trauma Overview
 - BioGears Trauma Models
- Medical Interventions
 - BioGears Treatment Models
 - What Can Go Wrong?
- BioGears Simulation Results
- Conclusions
 - How BioGears Can Impact This Field?



Multi-Trauma Physiology


- Multi-Trauma: subjected to multiple traumatic injuries
 - Hemorrhaging across body
 - Tension pneumothorax
 - Burn
 - Musculoskeletal injury
 - Traumatic brain injury
- Physiological changes vary not only based on injury, but on
 - Patient variability
 - Interventions
 - Time to intervention

Multi-Trauma Physiology




- Multi-Trauma in BioGears
 - Hemorrhage
 - Tension pneumothorax
- Tactical Combat Casualty Care (TCCC) emphasizes hemorrhage and airway management^[1]
- Goals
 - Medic: Decrease mortality and morbidity of the patient
 - BioGears: Physiological simulation to accurately depict multi-trauma scenarios

Incidences Background

- Hemorrhage: Bleeding, either internally or externally, from a broken blood vessel
 - May lead to inadequate blood gas exchange
 - Largest cause of combat deaths, over 80%^[1]
- Pneumothorax: When air gets trapped in the space between the lungs and chest wall (pleural cavity)
 - One way valve effect of respiratory exchange
 - Traumatic vs non traumatic
 - 10-15% of preventable deaths^[1]
 - third most potentially survivable cause of death^[2]
 - 5% of subjects arrived at a support hospital w/o a definitive airway^[1]




BioGears Trauma Models




- Hemorrhage
 - Circuit pathway with a switch connected to ground is modified to instantiate bleeding
 - Resistance based on flow rate
- Tension Pneumothorax
 - Closed – traumatic injury causing a “sucking chest wound”
 - Progressive build up of air within pleural space compartment of BioGears circuit structure

Input Diagram



Now What?

- If we can accurately depict a patient state after injury, we can also explore care/treatment scenarios



Background: Medical Intervention

- In prolonged field care, combat medics have limited tools
 - Substances
 - Tranexamic Acid (TXA)
 - Tylenol
 - Whole Blood (more commonly FWB is required)
 - Tools
 - Transfusion bag
 - Tourniquet
 - Decompression needle
- However, specificity and variability can make decision making hard
 - If medics aren't trained on specific physiology, patient response becomes harder to predict
- BioGears set-up reduces computational cost of modeling/simulation

Background

- **Tranexamic Acid (TXA)**

- Control of bleeding takes precedence over infusing fluids (TCCC) ^[2]
- Clinical Practice Guidelines^[2]
 - TXA should be given to casualties at risk of hemorrhagic shock as soon as feasible

- **Transfusion**


- Increase in blood volume increases blood gas transfer and immune health
- Massively transfused casualties have a high (33%) mortality rate^[2]
- In austere conditions, fresh whole blood is obtainable via a walking blood bank (WBB) program^[2]
 - Improperly done can lead to transfusion reaction

- **Needle Decompression**

- Decompress any suspected tension pneumothorax^[1,2]
- Casualties with multi-trauma showing no pulse or highly decreased respiration should have a bilateral needle decompression performed^[1,2]


Background: Transfusion Reactions

- Any adverse event associated with blood transfer
 - Hemolytic transfusion reaction (HTR) is most common [3]
- The blood type on identification tags is occasionally incorrect (last correlated data equated to about 4% inaccurate) [4]
- Can cause noticeable metabolic and respiratory distress in patients [3]



Tranexamic Acid: BioGears Validation

- Pharmacokinetics [5,6]




Tranexamic Acid: BioGears Validation

- Pharmacodynamics

- Short Scenario based on calculated blood volume after hemorrhage [7,8]
- Long Scenario based on time of survival after dosage [9]


1



2

Hemorrhage Survival Time	Hemorrhage w/ TXA Survival Time	Change	Time To Death
4691.4 s	5341.2 s	+ 13.9 %	Increases

BioGears Transfusion Model







Model: Transfusion Reaction


- Predator-Prey Relationship of Agglutination^[10]
- Assume:
 - Surface area estimated as cylinders
 - Four cells is “stable” and remove those cells/attachments
 - Only unlike cells can agglutinate

	Expected	Timing	Value	% Change	Result	Notes	Source
Core Temperature	Increase(>1 degC)	after	37.1	0.27027027		one degree change reached at 1820 secs	Davenport2005Pathophysiology/Anaesth2014Adverse
Respiration Rate	Increase		20.55	30.14566181			Davenport2005Pathophysiology
Blood Pressure	Change(Systolic, +/- ~20% @ severe)	during/shortly after	117	2.631578947		Can show increase OR decrease, key is change, severe not defined	Davenport2005Pathophysiology/Anaesth2014Adverse
Pulse	Increase(Tachycardia)	Immediately after	93.75	28.12628126		Tachycardia reached about 2.5 mins after end of transfusion	Davenport2005Pathophysiology/Anaesth2014Adverse

Agglutination Results (4 hrs)




Model: Needle Decompression




- BioGears incorporates a circuit component to permit airflow based on pressure differentials between pleural space and environment
 - Modified Resistance based on scaling factor
 - Appropriately tunes the circuit to vary breathing/air flow over time

Actions and Scenario Exploring




- Standard male patient
- Incidence
 - Closed left tension pneumothorax - moderate
 - Liver hemorrhage – Class III
 - Left arm hemorrhage – Class II
- Possible Interventions
 1. Do nothing
 2. Follow CPG
 3. Improper Care

Scenario Timelines







Standard Male Patient




Scenario 1: Do Nothing

- Exsanguination median time to death is 1.6 hours (5760 s) ^[12]






Scenario 2: Clinical Practice Guidelines



- — is traumatic incident: initiate multi-trauma
- — is hemostatic dressing and TXA injection: slowed blood loss
- — is needle decompression: restore breathing
- — is Blood transfusion to restore blood levels and relieve cardiovascular distress

Scenario 3: Medical Errors




- — is traumatic incident: initiate multi-trauma
- — is needle decompression: restore breathing
- — is TXA injection: low dose but still slowed blood loss
- — is Blood transfusion, stopped early due to HTR

Conclusions

- BioGears can accurately depict physiological responses to multi-trauma events and a variety of interventions
- Numerous applications
 - Training: Help less experienced personnel understand the dynamics of human physiology
 - Response to proper and improper care
 - Decision making: Using faster than real time simulations can help determine if an approach is more or less likely to help a patient

Acknowledgements

- A special Thank You goes out to:
 - The BioGears team for assistance in developmental stages of presented material
 - Austin Baird
 - Steven White
 - Matthew McDaniel
 - Lucas Marin
 - Jenn Carter
 - Jonathan Wingate, MD at the University of Washington for his review and expertise of the physiological components of presented material
 - Our research partners and mentors at TATRC and JPC-1
- BioGears would not be possible without their continued guidance, encouragement, leadership, and stewardship



References

1. Gerhardt, R. T., Mabry, R. L., De Lorenzo, R., & Butler, F. (2012). Fundamentals of combat casualty care. *Combat Casualty Care: Lessons Learned from OEF and OIF*. Washington, DC: United States Department of Defense, 85-120.
2. Butler, F. K., & Kotwal, R. S. (2017). Tactical combat casualty care. Version 5.
3. Strobel, E. (2008). Hemolytic transfusion reactions. *Transfusion Medicine and Hemotherapy*, 35(5), 346.
4. Cap, A. P., Beckett, A., Benov, A., Borgman, M., Chen, J., Corley, J. B., ... & Kane, S. F. (2018). Whole blood transfusion. *Military medicine*, 183(suppl_2), 44-51.
5. Pilbrant A, Schannong M, Vessman J. Pharmacokinetics and bioavailability of tranexamic acid. *Eur J Clin Pharmacol*. 1981;20(1):65-72. PubMed PMID: 7308275.
6. Nilsson, I. M. (1980). Clinical pharmacology of aminocaproic and tranexamic acids. *Journal of Clinical Pathology. Supplement (Royal College of Pathologists)*, 14, 41.
7. Kushioka, J., Yamashita, T., Okuda, S., Maeno, T., Matsumoto, T., Yamasaki, R., & Iwasaki, M. (2017). High-dose tranexamic acid reduces intraoperative and postoperative blood loss in posterior lumbar interbody fusion. *Journal of Neurosurgery: Spine*, 26(3), 363-367.
8. Gupta, K., Rastogi, B., Krishan, A., Gupta, A., Singh, V. P., & Agarwal, S. (2012). The prophylactic role of tranexamic acid to reduce blood loss during radical surgery: A prospective study. *Anesthesia, essays and researches*, 6(1), 70.
9. Roberts I, Shakur H, Coats T, Hunt B, Balogun E, Barnetson L, Cook L, Kawahara T, Perel P, Prieto-Merino D, Ramos M, Cairns J, Guerriero C. The CRASH-2 trial: a randomised controlled trial and economic evaluation of the effects of tranexamic acid on death, vascular occlusive events and transfusion requirement in bleeding trauma patients. *Health Technol Assess*. 2013 Mar;17(10):1-79. doi: 10.3310/hta17100. PubMed PMID: 23477634; PubMed Central PMCID: PMC4780956.
10. Dolgosheina, E. B., Karulin, A. Y., & Bobylev, A. V. (1992). A kinetic model of the agglutination process. *Mathematical biosciences*, 109(1), 1-10.
11. Davenport, R. D. (2005, July). Pathophysiology of hemolytic transfusion reactions. In *Seminars in hematology* (Vol. 42, No. 3, pp. 165-168). WB Saunders.
12. Cripps, M. W., Kutcher, M. E., Daley, A., McCreery, R. C., Greenberg, M. D., Cachola, L. M., ... & Cohen, M. J. (2013). Cause and timing of death in massively transfused trauma patients. *The journal of trauma and acute care surgery*, 75(2 0 2), S255.